

Lab Script

Pre labs

- Q6) Biological cell has internal compartments, radii $0.3-0.5 \mu\text{m}$. Estimate diffusion constant according to

$$\eta = 2 \times 10^{-3} \text{ kg / ms}$$

motion not hindered by other components.

$$\Rightarrow D = \frac{k_B T}{6\pi\eta R}, \text{ Let } T = 300\text{K}$$
$$k_B = 1.38 \times 10^{-23}$$

for $R = 0.3 \mu\text{m} = 0.3 \times 10^{-6} \text{ m}$

$$\hookrightarrow D = \frac{1.38 \times 10^{-23} \times 300}{6\pi \times 2 \times 10^{-3} \times 0.3 \times 10^{-6} \text{ m}} = \frac{4.14 \times 10^{-21}}{1.13 \times 10^{-8}} = 3.66 \times 10^{-13} \text{ m}^2/\text{s} = 0.366 \mu\text{m}^2/\text{s}$$

for $R = 0.5 \mu\text{m} = 0.5 \times 10^{-6} \text{ m}$

$$\hookrightarrow D = \frac{4.14 \times 10^{-21}}{6\pi \times 2 \times 10^{-3} \times 0.5 \times 10^{-6} \text{ m}} = 2.20 \times 10^{-13} \text{ m}^2/\text{s} = 0.220 \mu\text{m}^2/\text{s}$$

$$\therefore D \approx 0.22 \text{ to } 0.37 \mu\text{m}^2/\text{s}$$

- Q7) Reynolds Number Re for a cell with $L = 1 \mu\text{m}$ swimming $20 \mu\text{m/s}$ in a fluid with viscosity 10^{-3} Kg/ns & density 10^3 Kg/m^3 .
Cell is Turbulent or Laminar?

$$\Rightarrow Re = \frac{\rho L v}{\eta} \quad (\text{to find Turbulent or Laminar})$$

Given $L = 1 \mu\text{m} = 10^{-6} \text{ m}$, $v = 20 \mu\text{m/s} = 2 \times 10^{-5} \text{ m/s}$, $\eta = 10^{-3} \text{ Kg/ms}$
 $\rho = 10^3 \text{ Kg/m}^3$

$$\rightarrow Re = \frac{10^3 \times 10^{-6} \times 2 \times 10^{-5}}{10^{-3}} = \frac{2 \times 10^{-8}}{10^{-3}} = 2 \times 10^{-5} < 10^{-3}$$

(Low \Rightarrow Laminar)

Q8)

E coli sol. of optical density at 600nm of $OD_{600} = 1.0$
 This corresponds to 2×10^8 cells/ml!

Want ≈ 10 cells in FOV, dilute? By what?

\Rightarrow Expected number is $N = VC$

$$\begin{aligned} V_{\text{view}} &= A_{\text{FOV}} \times h \\ &= (N_x \times N_y) \times s \times h \rightarrow \text{Thickness} \\ &\quad \uparrow \quad \uparrow \\ &\quad \text{pixels} \quad \text{Calibration } (\mu\text{m}/\text{px}) \\ &= (1288 \times 964) \times 0.1 \mu\text{m}/\text{px} \times 5 \mu\text{m} \end{aligned}$$

$$\text{Field Size} \approx N_x s \times N_y s = 128.8 \mu\text{m} \times 96.4 \mu\text{m}$$

$$\text{Area} \approx 1.29 \times 10^4 \mu\text{m}^2 \times h$$

$$\text{Vol} = 6.21 \times 10^4 \mu\text{m}^3 = 6.21 \times 10^{-8} \text{ mL}$$

$$\therefore N = (2 \times 10^8) (6.21 \times 10^{-8}) \approx \underline{\underline{12.4}}$$

\therefore Yes dilute as now we see 12 cells

$$\text{dilution factor} \approx \frac{12.4}{10} = \underline{\underline{1.24}}$$

$$\therefore \underline{\underline{1.2 \times}} \quad \left(1 \text{ part culture \& } 0.2 \text{ part medium} \right)$$

1:2

Lab 2 Session 3: Microscopy and Cell Motility (Live Cells & Bacterial Motion)

Date: _____ Feb 2026 Lab Partner: Nathan Unhorn
Recorder: Ahilan Kumaresan

II. APPARATUS

Refer to Session 1, Section II (pg. 2)

Additional items for Session 3:

Item	Description
E. coli HCB1274	Tumble-disabled mutant strain
E. coli DH5 α (optional)	Wild-type strain
S. cerevisiae	Baker's yeast
HT1080 prepared slides	Fixed human cells (stained/ unstained)
Sugar, warm water	For yeast activation

III. VARIABLES

Type	Variable	Description
Independent	Cell type	E. coli, Yeast, HT1080
Dependent	$x(t)$, $y(t)$	Cell position over time
Dependent	v	Swimming velocity
Dependent	D_{eff}	Effective diffusion coeff
Control	Temperature	Room temp ($\sim 21^{\circ}\text{C}$)
Control	Frame rate	Match to cell speed
Control	Objective	100x oil immersion

LAB 1 Session 3

Lab 1 Session 3

I. GOALS

1. Observe and measure size of eukaryotic cells:
 - Yeast (*S. cerevisiae*): expect $\sim 5-10 \mu\text{m}$ diameter
 - HT1080 human cells: expect $\sim 20-50 \mu\text{m}$ diameter
2. Prepare *E. coli* (HCB1274 α mutant) sample with ~ 10 cells in FOV
3. Capture 60+ second video of bacterial motion at appropriate fps
4. Track bacterial trajectories using MTrack2
5. Compare bacterial motion to bead Brownian motion:
 - Is motion diffusive or ballistic?
 - Calculate effective diffusion coefficient
 - Compare MSD slope to beads
6. Determine if bacteria show directional (swimming) motion vs random Brownian motion

S. CELL SIZE - PROCEDURES and NOTES (Sec. 3.3, p. 9)

S.1 Yeast growth and wet-mount preparation

1) Prepared yeast culture in Falcon tube (refer to Sec. 3.3, p. 9):

Take clock and thermometer for room.

Tube volume: 15 ml Falcon

V_{water} (warm, body temperature) = 10 ml

m_{sugar} added = 2 g

m_{yeast} added = 2 g

Stirred gently for _____ s and left to sit for _____ min

Time when bubbling ~~STARTS~~ STABLES _____

(Write them like you did for Brownian motion; this is

Note: labscript footnote says amounts do not need to be exact, but I recorded what we actually used.

Cell size s

2) Wet mount (thin sample without spacer - "wet mount" in Sec. 3.3):

Independent variables:

Volume of yeast mixture on slide: _____ μ

Coverslip: _____ x 080) mm, #1 thickness (0.13-0.17 mm)

/ d

Any extra dilution? (yes/no)

If yes, yeast : water = 1 : _____

Dependent variables:

- Cell linear size along major axis

3) Brief notes on first visual check at low magnification (10x or 40x):

Density of cells: _____

Are individual cells resolvable _____

Any clumps / bubbles / debris? _____

5.2 Imaging yeast cells

Objective used for size measurement: _____x (write mag and NA if printed on lens)

Kohler illumination: set up as in Protocol:
Microscope Setup (same as previous session).

Image acquisition settings (Vision Assistant):

Camera: FLIR BlackFly BFS-U3-_____

Exposure time $t_c =$ _____ s

Gain or brightness settings: _____

Yeast image filenames

(store in Data/2026-__-__/Yeast/):

yeast_img_01.tif (path: _____)

Yeast

abcdefghijklmnopqrstuvwxyz

5.3 Yeast size measurement, formulas and slots

Calibration used here:

pixel size = _____ $\mu\text{m}/\text{pixel}$ \pm _____ $\mu\text{m}/\text{pixel}$

yeast_img_02.tif (path: _____)

(Copied from calibration section done previously with: _____
eyepiece magnification irrelevant, as emphasized in
labscript Sec. 3.1.)

Short notes per image (crow r

For each cell, measured number of pixels along major
axis:

Image 2 notes: _____

Image 3 notes: _____

Cell 1:

image = _____

N pixels major = _____ px

L major = N pixels major \times pixel size =
= _____ μm



Yeast

Cell 2

image = _____

N_pixels_major = _____

—
x

L_major = _____

—

m

C

I 3:

image = _____

N_pixels_major =

___ px

L_major = _____

S.4 HT1080 cell imaging (still Sec. 3.3, p. 9)

HT1080 slide description from labscript:

"HT1080 cells are epithelial cells derived from fibrosarcoma cells ... fixed mount, either unstained or stained blue."

Slide label and type used today:

Fixed HT1080, unstained / stained blue (circle one)

Label text on slide:

Imaging:

Objective: _____x

Camera and illumination settings similar to yeast (note any differences):

Exposure: _____s

Aperture setting: _____

Notes on contrast / color:

HT1080 images (store in Data/2026-__-__/HT1080/):

HT1080_img_01.tif (path: _____)

Shilaw
...

6. BACTERIAL MOTION - PROCEDURES and NOTES (Sec. 3.4, pp. 9-10)

6.1 Overview from labscript (write in notebook for context)

From Sec. 3.4 'Bacterial motion':

- Beads are passive (Brownian), bacteria are active materials with self-propulsion
- E. coli normally run at ~constant velocity for ~1 s, then tumble to choose a new direction
- Mutant strain HCB1274 has tumble disabled (run-only).
- Wild-type strain DHSα shows full run-tumble.
- "The experiment proceeds exactly as in Sec. B for polystyrene beads."

6.2 Dilution estimate (Prelab 8 logic)

Given (from labscript Prelab 8):

$$OD_{600} = 1.0 \quad 2 \times 10^8 \text{ cells/ml}$$

Goal: about 10 cells in field of view (FOV).

First, compute FOV volume:

$$N_x = 1288 \text{ px}, N_y = 964 \text{ px (camera)}$$

$P_{x \text{ Bac}}$ → pixel size Bacteria (Objective) = _____ $\mu\text{m/px} = \frac{\text{distance on micrometer (}\mu\text{m)}}{\text{no. of px}}$

Must Take multiple values & find uncertainty.

$$d_x = N_x \times P_{x \text{ Bac}} = 1288 \times \text{_____} = \text{_____} \mu\text{m}$$

$$d_y = N_y \times P_{x \text{ Bac}} = 964 \times \text{_____} = \text{_____} \mu\text{m}$$

$d_z \approx 5 \mu m$ (This is what Lab Script told us
Better ways to find This ???)

$$\begin{aligned}\Rightarrow \text{Vol. view} &= d_x \times d_y \times d_z \\ &= \text{---} \mu m^3 \\ &= \mu m^3 \times 10^{-12} = \text{---} \text{ mL}\end{aligned}$$

Now, for a target of ~10 cells in View

$$C_{\text{Need}} = \frac{N_{\text{Target}}}{\text{Vol. view}} \left(\frac{\text{cells}}{\text{mL}} \right)$$

$$\Rightarrow C_{\text{stock}} = 2 \times 10^8 \text{ cells/mL}$$

$$\text{dilution factor} = \frac{C_{\text{stock}}}{C_{\text{need}}}$$

oo We choose; $V_{\text{stock}} = \text{---} \mu L$, $V_{\text{Total}} = \text{---} \mu L$
our Dilution = $\frac{V_T}{V_{\text{stock}}} =$

Observation:

