

Lab 1 Session 6: Spatial Filtering (Fraunhofer Diffraction and Making Masks)

Date: 27 January 2026 Lab Partner: Nathan Unrhu

For Table of Contents, refer to session 4,
Pg 40

1. Session Goals

- i) Compare diffraction pattern to theory
- ii) Calibrate the Ronchi ruling (slit spacing a and slit width b)
- iii) Establish a quantitative baseline for spatial filtering 4. Complete spatial filtering mask cases (a), (b), (d), (e), (f), (j), (k)

2. Reference: Pages 39-43 (Part A: Fraunhofer Diffraction, Part B: Spatial

3. Variables

Independent: Mask configuration (orders transmitted) **Dependent:** Real-space image intensity profile, Fourier-space diffraction pattern

Control: Camera settings (gain, exposure, gamma), LED wavelength, grating period, optical alignment

4. Apparatus:

In Addition to the Material listed on Pg 2-3. We need

- i) Additional Translational Stage
- ii) Wires, Black Tape, Allen Key of different Size
- iii) A rectangular Frame with a Post and Saddle

5. Additional Background

6. Procedure:

A: Fraunhofer Diffraction

Following procedure on Lab Script pg. 39-40 3.1 Procedure 1.

- Replaced variable diffraction grating with 10 lp/mm Ronchi ruling microscope slide
- Set camera parameters: Exposure = 43000, Gamma = 1, Gain = 18
- Focused on "ZERO" lettering on Ronchi ruling for coarse focus
- Fine focused on blemishes (chipped metal on lettering) for optimal sharpness
- Inserted 10 nm bandpass filter in front of LED
- Verified Fourier image centro-symmetry
- Recorded real-space and Fourier-space images
- Extracted line profiles for calibration measurements

Observation:

Despite turning the dial on the condenser lens, the focus appeared to change unexpectedly. Fine focusing required locating small defects (chipped metal) on the ruling surface rather than the grating lines themselves. This provides sharper focus targets than the periodic structure due to the Talbot effect creating multiple focal planes for periodic objects.

Observation:

The Fourier image was centro-symmetrical, confirming that the optical elements are properly aligned square to the optical axis. This symmetry verification is essential before quantitative measurements.

Calibrate the Real image:

10 lines on the camera is from 602px \rightarrow 1150px

Pixel span 1150 - 602 = 548 pixel for 10 period

$$\text{Scale factor} = \frac{a}{\text{Pixel per period}} = \frac{100 \mu\text{m}}{54.8 \text{ px}} = 1.82 \mu\text{m} / \text{pixel}$$


With the Real space camera having a pixel size $\Delta_{\text{pixel}} = 3.45 \mu\text{m}$.

$$M_{\text{real}} = \frac{\text{image scale}}{\text{pixel size}} = \frac{54.8 \times 3.45 \mu\text{m}}{100 \mu\text{m}} = \frac{189 \mu\text{m}}{100 \mu\text{m}} = 1.89 \times$$

$$\text{Expected magnification} = \frac{f_5}{f_4} = \frac{300}{150} = 2.0 \times$$

The values are close, but the difference is from the lens position not being exactly at focal length

Calibration for Fourier Image:



Your label	Pixel distance	Ratio to 1st	Actual order m	Value
"1st"	78 px	1.00	$m = 1$	200 (1600, Expected, 1st, 19 Grain) 100
"2nd"	238 px	3.05	$m = 3$	
"3rd"	396 px	5.08	$m = 5$	
"4th"	553 px	7.09	$m = 7$	

Value

800 Expected 1200 ps
120
50
100
50

We Note that we only see odd orders $m = 1, 3, 5, 7$.

This is what we observe from a 50% duty Cycle Ronchi ruling.

The even orders are suppressed because $m = 2, 4, 6$, because $\sin\left(\frac{m\pi b}{a}\right) = \sin\left(\frac{m\pi}{2}\right) = 0$ for even m when $\frac{b}{a} = \frac{1}{2}$

Feynman Camera Calibration

$$\theta_m = \frac{m\lambda}{a} = \frac{m \times 525 \times 10^{-9}}{100 \times 10^{-6}} = m \times 5.2 \text{ mrad}$$

$$K_{cal} = \frac{\theta_1}{\Delta p_1} = \frac{5.25 \text{ mrad}}{78 \text{ px}} = 0.0673 \text{ mrad/px}$$

Verification: Measured vs Theory			
Order m	Measured (px)	Theory: $m \times 78 \text{ px}$	Discrepancy
1	78	78	0% (reference)
3	238	234	+1.7%
5	396	390	+1.5%
7	553	546	+1.3%

Excellent linearity! The small systematic offset ($\sim 1.5\%$) could be due to:

- Slight deviation from small-angle approximation at higher orders
- Minor lens aberrations
- Grating period slightly different from nominal $100 \mu\text{m}$

\Rightarrow Extracting Grating Period

$$\Delta p = \frac{78 + 238/3 + 396/5 + 553/7}{4} = \frac{78 + 79.3 + 79.2 + 79.0}{4} = 78.9 \text{ px/order}$$

$$a = \frac{\lambda \cdot \text{fobj}}{\Delta x} = \frac{\lambda}{\theta_1} = \frac{525 \text{ nm}}{5.25 \text{ mrad}} = 100 \mu\text{m}$$

Summary: Both Calibrations

Method	Parameter	Value
Real-space	Scale	1.82 $\mu\text{m}/\text{pixel}$
	Magnification	1.89 \times
	Period a	100 μm (known)
Fourier-space	Calibration	0.0673 mrad/pixel
	Period a (verified)	100 μm ✓
	Duty cycle b/a	~ 0.5 (even orders absent)

\Rightarrow Pixel Value


B. Spatial Filtering


Following procedure on Lab Script pg.
41-43

2:57 PM

- 1. Verified Fourier patterns are symmetrical and well focused
- 2. Placed mask post at Fourier plane (after condenser lens, before beam splitter)
- 3. Mounted translational stage for precise mask positioning

Saved as file: Default and blocked-ve1to1







At 1

At the Fourier Plane, after the condenser lens, before the Beam Splitter. We placed a post, on this post we are placing masks to do diffraction

TIP: We are using a hellend so translation works very very well.



Placing the focus mask after the condenser lens. at the fourier plane

(a) All orders (reference case)
Saved Images as "RonchiReal-A.tiff"
"RonchiFourier-A.tiff"
Default image. Looks bright and the real image is all well illuminated

We accidentally increased the Irirs so new images as


Saved Images as "RonchiReal-A2.tiff"
"RonchiFourier-A2.tiff"

(b) Zero order only;
Intensity Reduced



(b) Zero order only;
Intensity Reduced

Saved Images as "RonchiReal-B.tiff"
"RonchiFourier-B.tiff"



(c) +1 only;
(d) 0, +1 (Abbe minimal criterion);

We have moved it forwards, to make the mask larger and slid it on the translational stage.

Saved Images as "RonchiReal-D.tiff"
~~Do not~~ "RonchiFourier-D.tiff"

Following procedure in B-D
More multiple lines in-between each solid white line

m_{\max}	Orders Transmitted	Iris Setting
1	-1, 0, +1	Smallest (just fits ± 1)
3	-3, -1, 0, +1, +3	Medium (just fits ± 3)
5	-5 to +5	Larger (just fits ± 5)
7	-7 to +7	Even larger (just fits ± 7)
All	Reference	Fully open

(f) +1, -1;

Using the thinnest Allen key and a similar size Screw driver.

Saved Images as "RonchiReal-F.tiff"
"RonchiFourier-F.tiff"

(g) +2, 0, -2 and also +2, -2 (note: for a perfect square wave, even harmonics such as ± 2 would be suppressed);

- (h) All positive orders (oblique dark field);
- (i) All positive orders and the zero order (half shadow);
- (j) All orders except zero (dark-field / dark-ground method);

Using a Slit Mask, just blocked the Central Maximum
Saved Images as "RonchiReal-F.tiff"
"RonchiFourier-F.tiff"

Just a uniform dark screen on the Real image.

(k) $-m_{\max}, \dots, -1, 0, +1, \dots, +m_{\max}$
(low-pass filtering; choose at least two values of m_{\max} in

addition to 1 e.g. 1, 3, 5);

Doing outer Filtering, from -2 to + 2.

And

Saved Images as "RonchiReal-K2.tiff"
"RonchiFourier-K2.tiff"