

Lab 1 Session 4 SF

Prelab Question 3.

1. Let the grating be periodic in x with period d . Over one period, model its amplitude transmittance as a rectangular pulse:

$$t(x) = \begin{cases} 1 & , 0 \leq x \leq \omega \\ 0 & \omega \leq x \leq d \end{cases}, \text{ periodic with period } d$$

$$\Rightarrow \text{Duty Cycle; } f = \frac{\omega}{d} \quad (0 < f < 1)$$

Fourier Series form:

$$t(x) = \sum_{n=-\infty}^{\infty} C_n e^{i 2\pi n x / d}$$

$$\begin{aligned} \text{& fourier coeff: } C_n &= \frac{1}{d} \int_0^d t(x) e^{-i 2\pi n x / d} dx \\ &= \frac{1}{d} \int_0^{\omega} e^{-i 2\pi n x / d} dx \end{aligned}$$

$$\text{DC offset term, } C_0 = f$$

for $n \neq 0$

$$C_n = \frac{1}{d} \left[\frac{e^{-i 2\pi n \omega / d}}{-i 2\pi n / d} \right]_0^{\omega} = \frac{-1}{i 2\pi n} \left(e^{-i 2\pi n f} - 1 \right)$$

$$C_n = \frac{\sin(\pi n f)}{\pi n} e^{i \pi n f} \quad (n \neq 0)$$

- 2) In a 4f imaging system, fourier plane contains discrete diffraction orders n with complex amplitudes proportional to C_n . A spatial filter modifies these orders by a complex factor M_n (amplitude & phase)

Filter is modelled as:

$$E(x) = \sum_{n=-N}^N (C_n M_n) e^{j 2\pi n x / d}$$

Absorbal, $I = |E(x)|^2 = E(x) E^*(x)$
 \approx complex conjugate

Code Below:

Pre-lab Question 3

January 20, 2026

Ahilan Kumaresan

```
[4]: import numpy as np
import matplotlib.pyplot as plt

[15]: def ronchi_coeffs(N, duty_cycle):

    f = duty_cycle
    n = np.arange(-N, N+1, dtype=int)

    c = np.zeros_like(n, dtype=complex)

    # n=0
    c[n == 0] = f

    # n != 0
    nz = (n != 0)
    nn = n[nz].astype(float)
    c[nz] = (np.sin(np.pi * nn * f) / (np.pi * nn)) * np.exp(-1j * np.pi * nn * ↵f)

    return n, c

def make_mask(n, mode="", **kwargs):

    M = np.ones_like(n, dtype=complex)

    if mode == "":
        return M

    if mode == "block_dc":
        M[n == 0] = 0
        return M

    if mode == "keep_orders":

        keep = set(kwargs.get("keep", []))
        M[:] = 0
```

```

    for k in keep:
        M[n == k] = 1
    return M

if mode == "lowpass":
    nmax = int(kwargs.get("nmax", 1))
    M[np.abs(n) > nmax] = 0
    return M

if mode == "attenuate":
    scale = kwargs.get("scale", {})
    for k, s in scale.items():
        M[n == int(k)] *= complex(s)
    return M

if mode == "phase_shift":
    phase = kwargs.get("phase", {})
    for k, ph in phase.items():
        M[n == int(k)] *= np.exp(1j * float(ph))
    return M

raise ValueError(f"Unknown mode: {mode}")

def field_and_intensity(d, duty_cycle, N, x, mask=None):
    n, c = ronchi_coeffs(N, duty_cycle)

    if mask is None:
        M = np.ones_like(c, dtype=complex)
    else:
        M = mask

    # Build exp(i 2 pi x/d) efficiently: shape (num_orders, num_x)
    phase = np.exp(1j * 2*np.pi * np.outer(n, x) / d)
    E = np.sum((c * M)[:, None] * phase, axis=0)
    I = np.abs(E)**2
    return n, c, M, E, I

def normalize(I):
    return I / np.max(I)

# PARAMETERS
d = 100e-6 # period [m] (example: 10 lp/mm => d=0.1 mm = 100 μm)
f = 0.50      # duty cycle (open fraction); use your measured/assumed value
N = 15        # truncation order
x = np.linspace(-2*d, 2*d, 4000) # evaluate over several periods

```

```

# EXAMPLE
cases = [
    ("No filter",                      make_mask(np.arange(-N, N+1), mode="")),
    ("Low-pass abs(n)<=1",            make_mask(np.arange(-N, N+1), mode="lowpass", nmax=1)),
    ("Block DC (dark-field-ish)",      make_mask(np.arange(-N, N+1), mode="block_dc")),
    ("Keep only 1 (edge-ish)",         make_mask(np.arange(-N, N+1), mode="keep_orders", keep=[-1, 1])),
    ("Atten + / - 1 to 0.2",          make_mask(np.arange(-N, N+1), mode="attenuate", scale={-1:0.2, 1:0.2})),
]

```

A4_WIDTH = 8.27
A4_HEIGHT = 11.69

```

fig, axs = plt.subplots(4, 3, figsize=(A4_WIDTH, A4_HEIGHT),
                      layout='constrained', sharex=True, sharey=True)

for (label, M), ax in zip(cases, axs.flat):
    n, c, M, E, I = field_and_intensity(d=d, duty_cycle=f, N=N, x=x, mask=M)
    ax.plot(x*1e6, normalize(I))

    ax.set_title(label, fontsize=10)
    ax.set_xlabel("x [μm]", fontsize=8)
    ax.set_ylabel("I/max(I)", fontsize=8)
    ax.grid(True, alpha=0.3)

for ax in axs.flat[len(cases):]:
    ax.set_visible(False) # Hiding other non printing plots

fig.suptitle(f"Optical Fourier Filtering (N={N}, f={f})", fontsize=14)

plt.show()

```




Figure i-v)

Lab 1 Session 4: SF (Completing Setup 1 & Testing Image)
Date: 20 - Jan - 2026 Lab Partner: Absent

Table of Contents (For Session 4 & 5)

I. Session 4 (completing Setup 2 and Testing 1 and 2)

1. Goals
2. Apparatus
3. Sketch
4. References
5. Useful Background
6. Procedure
7. Testing
8. Post lab reflections and sources of errors.

II. Session 5 (Redo Test Two and Refine Set up for Data Collection)

Table of Contents:

1. Goals
2. Background/Theory
3. Procedure:
4. Post lab analysis

1. Goals

- Capture an image for the real panel, also one image of the real panel after the beam splitter is installed. For seeing the reduced intensity.
- Adjust exposure to have better contrast and visibility in the real plane.
- Complete Lab Set up 1 and Testing

2. Apparatus

Optics / mounts

- Beam splitter
- 90° short optical rail + feet (leveled)
- 2 irises/apertures for alignment on 90° rail
- Fourier imaging lens: 50 mm, $f = 100$ mm
- Lens tube for Fourier camera
- (Later) transverse translation stage + iris for Fourier plane filtering (confirm clearance)

Cameras

- Fourier camera (FLIR Blackfly): mounted on transverse translation stage (XY)
- Real-space camera already aligned from Setup Part 1


Tools

- Alignment cap with central dot
- Viewing card (white card / paper)
- Hex drivers, ball driver, post collars, spare saddles/posts (3", 4", 6")

Software

- LabVIEW: Dual_Camera_Line_Profile_v2.vi
- NI Vision Assistant (optional, but close it before LabVIEW)

3. Sketch:



Sketch 1: Top View of Secondary Rail & Elements Setup

4. References:

- Optical Systems Guide J.M. McGuirk
- An Introduction to Optics. Lipson
- Spatial Filtering Lab Script

5. Useful Background

Field Plane Conjugates (image of light source structure):

Field Stop Object Real-space Camera

Field stop controls illuminated area (field of view)

Changing field iris changes FOV but NOT resolution

Aperture Plane Conjugates (image of source angular distribution):

Aperture Stop Fourier Plane (determines resolution)

6. Detailed Procedure

i) Verified Beam Splitter Insertion and 90° Rail Alignment

Place beam splitter in main beam

- Verified beam splitter is approx. 22 cm downstream of objective lens (L4).

Verify:

This spacing must leave room between BS and Fourier plane for a transverse translation stage with iris (spatial filter later).

- visually confirm clearance: stage can mount on a saddle and sit near BS without collision.

Record

- Objective lens reference point: $104 \pm 0.05 \text{ cm}$
- BS placed at: $82 \pm 0.05 \text{ cm}$ downstream of objective (target 22 cm)
- Clearance check for future iris/stage: /

Notes:

There is space but I do not see the Fourier Plane image

Issue:

ii) Rotate BS (using a collar) so that the reflected beam propagates straight on one arm and away from the edge on the other.

- Place the smaller rail with the secondary camera perpendicular to the beam splitter, so that the ray from the beam splitter hits the camera.
- Use Iris and Aperture to verify the centre.

USE the white card and the raw from the BS to see and align the Camera Sensor to the light.

Expected:

I expected to see a diffraction pattern, just as before in page ()

Saved Image at Session4/Notes/

Name: Fourier-Image-Pinhole.tiff




Image 1: Fourier Image Through a pin hole.

Could not see a diffraction pattern, maybe we are out of focus? There seems to be some blurring. Perhaps the BS is unclean, too.

However, Success! As I slide the card along the rail, the beam always stays within the Secondary Rail (with secondary camera)

1. Fourier imaging lens and Fourier camera
2. Placed the Fourier Lens (FL) (50 mm diameter, $f = 100$ mm) in the secondary rail. Ensure correct placement by facing the side with the infinity symbol towards the BS.

Check
2:20

This was the last lens available, so it matches the description, giving confidence that our lenses are being used correctly.

Check :

iii) Verified that the fourier camera is on the translational stage and added a Lense tube in front to block stray light from the room.


Imagine the Fourier plane onto the Fourier camera so that up to the 7th diffraction orders fit on the CCD for a 10 Lp/mm grating on the test object. ()

Adjusted the Fourier Lens (FL) along the rail until we finally saw the diffraction pattern, for Best Focus!!
Yay!! FL pos: 29cm along secondary rail

Manually adjusted the railing a little, just the centre, to diffuse the pattern better.

Saved Image at Session 4/Notes/
Name: diffraction-pattern-secondary-rail-01.tiff.

Image 2: Diffraction Pattern from secondary rail Camera



Issue: Can only see 4 diffraction orders, There is probably more as when we adjust the mirror horizontally, we can see more but. How do I make all of them visible?

I tried to full the Camera and Lense Set up towards the BS, that did not work too well, it only got zoomed in.


Asked Zack (TA) for help, the idea was to either pull as we did before or take multiple images with different exposure times to get the lower intensity images.

something new was the error happening in adjusting the Exposure,

error message: "Error - Occured at Property Node"

Screenshot saved at Session4/Notes

Name: error-changing-exposure.png



3:00

The fix is: go to BOTH camera Acusiation in the Block Diagram for the Preset and switch off the Auto Gain and Exposure. We had done it.

Success: The diffraction patter is Centro Symmetrical, Sharp and Can see multiple orders.

iv) Final Notes for Set up 2:

Goal Met: We have a dual-camera spatial filtering setup. LabView shows the Fourier and Image Plane.




Image 4, Our Optics Setup

Added cross sections on both the Panels and explored the Live intensity profile.

Save images and a screenshot of the setup.

in Session 4/Notes

Name: labView-Set up-Screenshot.png

That is the line
Intensity Pattern!
Wow!




Image 5, Our LabView Setup

7. Testing

i) Test 1

Observe 3 things simultaneously

1. The four Camera images
2. The real-space image
3. With a card, the illuminated region or object

3:25 PM

Did test 1 on (pg 36) of Lab Notebook.


⇒ Varying field stop (field iris), pupil size

Expectation

⇒ We observed area on object & camera field of view change on images.

⇒ On Fourier Camera: Spot Size stays same but Intensity increases with more light let in.
We actually see more orders & more blurring actually

image: Field Iris



Saved images in Session 4/Notes

Saved 3 images of 3 different heuristic pupil size.

"Smallest-pupil-field-iris.tiff", "1-10th-pupil-field-iris.tiff"
& "2-10th-pupil-field-iris.tiff"

New Approach to Note Taking:
Zach, our TA recommended only referring procedure from our sources (like Lab Script) and writing only that is new and our personal notes!
This changes everything!
Ahilan Kumaresan

ii) Test 2

* Test 2 procedure on (pg 36) of lab Script

Varying Aperture Stop using an iris at Aperture stop plane.




The Zeroth order on the Four Camera

is at pixel (713,508)

↳ Marked the Zeroth order with a line

↳ Using a adjoint Saddle & remove the
(AS) pinhole & insert an iris at aperture stop plane

As usual make sure it is square & wing a up
stream pin hole; it is centred.



Aperture Stop (Pin hole)
Image 7:

4:06 PM

Facing an **error** because the new iris is not fully on the centre of the screen.

Looking at the light train, looks like the problem might be down stream too. Using a white card & another pinhole.

Trying to fix light train, making it aligned to optical rail.

Put Back the pinhole & trouble shooting.

Resolved:

4:21 PM

Finally!! The issue was on (L³) Condensor lens

Adjusting Secondary Rail elements for centred lens.
new centre (709, 620) pixel.

How ever something changed, its more teardrop shaped
Sand image in Session 4/ Notes


Name : new-diffraction-pattern.tiff (See below)

8. Post lab Reflection :

I made some significant headway, a lot of the accumulated errors were resolved :

- 1) Optics element mis alignment
- 2) Errors in Lab View Exposure Settings

However Some new errors popped up :



old diffraction patte

New diffraction patter

Probably this shows some unclean lenses?

i) Goal Standing

Goal 1: Capture real-space image before/after beam splitter MET ✓

- Successfully observed intensity reduction after BS insertion
- Adjusted exposure to compensate

Goal 2: Adjusted Exposure and got a better understanding of exposure settings. MET

Goal 3: Take image of real plane and upload it.
Yes Done.

Goal 4: Complete Lab set up 2 and Testing.
Partially. Testing to be redone next lab.

ii) Learnings:

1. Alignment is critical - small misalignments cause visible artifacts
2. The Fourier plane location can be found by sliding a card along the beam path - first you see real-space image, then diffraction patt
3. LabVIEW errors often require checking BOTH camera setting
4. The intensity cross-section tool is very useful for quick alignm

iii) What I Would Do Differently:

Record more quantitative measurements (pixel positions, peak widths)

- Take systematic images at each alignment step for comparison
- Note iris/aperture settings when taking images

iv) Questions to Investigate in Session 5:

1. Why are diffraction spots teardrop-shaped instead of circular?
2. Can I see all 7 orders for 10 lp/mm as expected?
3. What happens when I vary the field iris vs aperture iris?

Lab 1 Session 5:

Spatial Filtering (Finishing Setup 2 & Improving diffraction image resolution)

Date: 22-Jan-2026

Lab Partner: Nathan Unku

FOR TABLE OF CONTENTS, REFER TO SESSION 4 pg 40

1. Goals:

- i) Labscript procedure testing after setup
- ii) Error correction during setup
- iii) Understand how the Fourier plane spots are formed and how "sharpness" or "detail" shows up as one spot while brightness shows up as another

2. Background Information Physical intuition (no math)

- i) Light behaves like a wave. A wave cannot make an abrupt change without involving many different directions of motion.

11:46 AM

- ii) A sharp edge is like telling the wave: "Go from zero to full amplitude instantly." The only way a wave can do that is by combining many plane waves at different angles.

- Few angles smooth variation
- Many angles sharp transitions;

iii) Slit intuition (classic example)

- A wide slit ;(smooth spatial variation) narrow diffraction pattern
- A narrow slit ;(sharp spatial confinement) wide diffraction pattern

This is the same reason:

- Fine details large angles
- Coarse features small angles

iv) Filtering methods:

- Low pass filter: removes noise
- High pass filter: sharpness improvement
- Band-pass filtering: specific frequencies selected to show texture/patterns otherwise hidden by bright or sharp detail
- Directional filtering: orthogonal structures removed (only vertical or horizontal features)
- Notch filtering: removing specific noise
- Phase-contrast filtering: shifting frequencies to make transparent objects more visible
- Dark-field filtering: Reduce white highlights to see dark/small features
- Amplitude filtering: changes contrast
- Spatial frequency weighting: custom contrast for spec. freq.
- Order selection: resolution control
- Spacial differentiation: for edge/feature detection
- Holographic and correlation filtering (advanced): pattern recognition in real time

3. Detailed Procedure

i) Image Doubling: Trouble Shooting

We got a suggestion to use a filter paper that is translucent. Go through the different elements and try to determine which element starts the double first.

The issue could be the back reflection from one of the shiny parts of the elements.

However, My hunch is some misalignments of the equipments.

Looking at the current state, maybe this issue isn't too bad. The intensity pattern definitely show multiple values. But Maybe this isn't too bad ???


Saved as initial-state - diffraction pattern.tif

Observation

I can say that because there is both more dots on the fourier image but also, the intensity patterns send multiple smaller peaks.

These are not diffraction orders.

Image 1: Initial State Screenshot



Clearly There are multiple reflections

ii) Correction of Lense

1:42PM

Correction of error (image is off center and doubling of forrier plane spots)

- Starting from test step 2 in setup 2 of labscrip (vary the aperature)
- Placing a paper filter behind condensor lens fixed the doubling, concluding that the doubling was cause by misalignment of either incoming light onto condensor or condesor alignment. Realigning the condensor lens should correct this.
- Fixed doubling by pitching lens forward on track and right very slightly.




Figure 1: Ray sketch of Back Reflection

iii) Correcting Fersnel Fridnges

New Issue: Fersnel Fridnges on the real Image, This makes it clear that our has been image is out of focus! This could explain the doubling of our diffraction pattern! So my previous hunch of an "Unclean lense" is incorrect.

Saved this image as: fersnel-frindge-issue.tiff under Session5/Notes

Observation

Noticed how there are fringes around the letters we are focusing on.

Previously we thought this was insignificant but this seems to matter more for resolution as we move towards more tests.




Image 2: Fresnel Fringes on Image

Tip

When trying to focus & eliminate fresnel fringes try to focus on a small defect on the test object & make it sharp.

iii) Continuing with Test 2 again

- added the aperture stop before the forrier plane splitter, closing the aperture blocked the outer diffraction orders which blurs the image (outer orders at higher frequency contain the sharp detail information)
- Position of 0-order: 819, 655 pixels (zeroth-order-image.tiff) saved at session 4 notes
- Removed the pinhole temporarily and added the iris aperture in its place. This changed the brightness of the image without changing the illumination area much.
- Moved the small iris back toward the LED source and make the image edges sharper
- Saved Image: Correct-iris-position.tiff

Image3:

• Corrected sharp real
Image after moving iris

Compared to our
previous image the
edges are way
sharper!




Image4: Showing an example
of a defect on grating we
focused on



- inserted 10nm band filter in front of LED, eliminated the non-zero order diffraction patterns (collapsed to single frequency peak with black image)

iv) Resolution Test:

started with 10 lp/mm grading. With increased exposure and 10nm filter image is clearer (original wavelength is 525 +35nm)

Recorded images of 12 and 16.67lp/mm. Spacing between 16 (0th \rightarrow 1st order) is 202 pixels. Spacing for 12 is 132 pixels

(space below for calculation of magnification level)

$$\theta_1 = \sin^{-1}\left(\frac{\lambda}{f}\right) = \left(\frac{525 \cdot 10^{-9}}{16.67 \cdot 1000}\right)$$

$$= 0.501^\circ$$

$$202 \text{ pxl} \div 0.501^\circ = 403 \text{ pxl/deg.}$$

$$f_{\text{eff}} = \frac{\Delta x}{\theta_1} = \frac{202 \text{ pxl} \cdot 3.45 \mu\text{m}/\text{pxl}}{0.501^\circ}$$

$$= 139 \text{ mm}$$

v) Resolution test

- chose 26 and 50 lp/mm (easy and difficult to resolve)
- Prediction: closing the apperature iris slightly should block the higher order diffraction spots tied to the finner line spacing, blurring that section without blurring the coarse line spacing.
- Added an iris at the forrier plane. Reducing aperature eliminates higher order diffraction spots, reducing fine detail first (high frequencies = finest detail = farthest diffraction orders)
- Closed the iris halfway to make the finner grading lines blurred while the coarser ones still sharp (insert image of blured fine lines and sharp coarse lines)

Second order ↑ 1st order ↓ Zeroth order ↑ First order ↓

Image 5 :
26 lp/mm





Image 6 :
50 lp/mm



Image 7 :
26 and 50 lp/mm



All files saved under Session 5/Notes

Name : "26lp mm - fourier.tif"
"50lp mm - fourier.tif"
"26-50lpmm - fourier.tif"

Observation: We notice how in the 26lp/mm Grating, the first Order is very clear and bright. While Higher Orders are dim

On the 50 lp/mm we observe the second order bright while the first order is relatively dim.

We can see that when we position the grating such that the feild of view has 26lp/ mm on one side and 50 lp/mm on the other, we get a combined patter! This is pretty cool!

Success! :

confirmed overlap between fourier images of 26 and 50 lp/mm diffraction orders.

5) Post Lab Reflection and Suggestions for next Lab

The Diffraction Pattern (Page 63) was illuminating.

- 1) We were not able to see higher diffraction orders.

A solution we recommend is to tilt the Fourier camera so that the 0th order goes to one end. That way, we could imagine more orders.




Figure 2: Fourier Image Idea

- 2) While our image light illumination has dark spots, it is decently in focus after making some alignment corrections, and we are now able to see the effects of altering the fourier diffraction orders on the image results. We achieved our intial goals of correcting the focus and alignment of the fourier plane peaks. Our setup is sufficient to proceed with Ronchi ruling next session.

- 3) This lab fundamentally changed how I conceptualize image formation. Rather than thinking of lenses as simply "focusing light," I now understand them as performing spatial Fourier transforms, with the image plane representing the inverse transform of filtered frequency components. The ability to physically manipulate frequency content by placing objects in the Fourier plane --- and immediately see the effect on the image --- made abstract mathematical concepts concrete and intuitive.

Lab Notebook Submission Summaries

Spatial Filtering Lab (Lab 1) – Sessions 4 & 5

PHYS 332

Week 1 Submission: Session 4

Session Overview

Date: 20-Jan-2026 **Lab Partner:** Absent **Pages:** 40–53 **Objective:** Complete Setup 2 by adding beam splitter and secondary (Fourier) camera arm; perform initial testing of field and aperture stop effects.

[View Next Page: Format Error](#)

Item	Page
<i>Preparation</i>	
Table of Contents for Sessions 4 & 5	40
Goals: capture real-space images before/after BS; adjust exposure; complete Setup 2	41
Apparatus List (optics, cameras, tools, software)	42
Experimental Setup Sketch (top view of secondary rail configuration)	43
References and Background: Field & Aperture Plane Conjugates	43–44
<i>Procedure & Data Collection</i>	
Beam Splitter Insertion & Rail Alignment — placed BS at 82 ± 0.05 cm	44
Secondary Rail Setup — rotated BS, used white card to trace beam path	45
Fourier Lens Placement ($f = 100$ mm) — adjusted until diffraction pattern focused	45–46
Diffraction Pattern Acquisition — initially only 4 orders visible, adjusted mirror	46–47
LabVIEW Error Resolution — disabled Auto Gain/Exposure in Block Diagram	47
Final Setup Verification — dual-camera system confirmed operational (Images 4, 5)	48
<i>Testing</i>	
Test 1: Field Stop Variation — varied field iris, observed FOV change on real-space	49
<i>camera while Fourier spot size stayed constant (intensity changed)</i>	
Test 2: Aperture Stop Variation — zeroth order at pixel (713, 508)	50
Troubleshooting: Iris centering issue traced to condenser lens (L3) misalignment	50–51
<i>Post-Lab Analysis</i>	
New Issue: Diffraction spots became teardrop-shaped after adjustments	51
Errors Resolved: (1) optical element misalignment, (2) LabVIEW exposure settings	51–52
Goal Status: Goals 1–3 MET; Goal 4 (testing) PARTIAL	52
Key Learnings (alignment, Fourier plane location, LabVIEW, intensity cross-section)	52
What I Would Do Differently; Questions for Session 5	53

Key Files Saved:

- Fourier-Image-Pinhole.tiff
- diffraction-pattern-secondary-rail-01.tiff
- error-changing-exposure.png
- LabView-Setup-Screenshot.png
- new-diffraction-pattern.tiff
- Smallest-pupil-field-iris.tiff
- 1-10th-pupil-field-iris.tiff
- 2-10th-pupil-field-iris.tiff

Week 2 Submission: Session 5

Session Overview

Date: 22-Jan-2026 **Lab Partner:** Nathan Unhm **Pages:** 54–64

Objective: Resolve image doubling and focus issues from Session 4; understand physical relationship between Fourier plane spots and image sharpness; perform resolution tests with multiple gratings.

Item	Page
<i>Preparation & Background</i>	
Goals: complete testing, correct errors, understand sharpness vs. brightness in Fourier plane	54
Physical Intuition: wave behavior, sharp edges require many plane wave components	54–55
Slit Diffraction: wide slit → narrow pattern; fine details → large angles	55
Filtering Methods Overview (11 methods documented)	56
<i>Procedure & Troubleshooting</i>	
Image Doubling Diagnosis — used filter paper to trace source; multiple intensity peaks	57–58
Resolution: misalignment of light onto condenser; fixed by pitching lens forward/right	58
Fresnel Fringes Issue — fringes around text indicate out-of-focus image	58–59
Focus Tip: focus on small defect on test object and make it sharp	59
<i>Testing & Data Collection</i>	
Test 2 Continued: Aperture stop blocks outer orders → blurs image (high-freq info lost)	60
10 nm Band Filter Insertion — eliminated diffraction spreading, cleaner pattern	61
Resolution Test (10, 12, 16.67 lp/mm) — measured pixel spacings between orders	61
Diffraction Angle Calculation: $\theta_1 = \sin^{-1}(\lambda/d) = 0.501 = 0.00874$ rad	61
Resolution Test (26, 50 lp/mm) — predicted and confirmed selective blurring	62–63
<i>Closing iris halfway: fine lines (50 lp/mm) blur, coarse lines (26 lp/mm) stay sharp</i>	
Combined Grating Pattern — positioned FOV to capture both frequencies simultaneously	63
<i>Post-Lab Analysis</i>	
Reflection: diffraction pattern images illuminating; could not see all higher orders	64
Suggestion: tilt Fourier camera so 0th order at edge to fit more orders on sensor	64
Key Files Saved:	

- initial-state-diffraction-pattern.tiff
- fersnel-fringe-issue.tiff
- Correct-iris-position.tiff
- zeroth-order-image.tiff
- 26lp-mm-fourier.tiff
- 50lp-mm-fourier.tiff
- 26-50lpmm-fourier.tiff

Key Measurements:

Parameter	Value	Notes
LED wavelength	525 nm	Bandwidth 35 nm
Zeroth order (initial)	(819, 655) px	Before correction
Zeroth order (corrected)	(709, 620) px	After condenser fix
16.67 lp/mm spacing	202 pixels	0th to 1st order
12 lp/mm spacing	132 pixels	0th to 1st order
θ_1 (16.67 lp/mm)	0.501°	Calculated

Equipment: Fourier lens (50 mm dia, $f = 100$ mm), FLIR Blackfly camera on XY stage, 10 nm band filter, resolution target (10–50 lp/mm gratings).