

Prelab Lab 2

February 3, 2026 3:33 PM

If the microscope's calibration factor is 100 nm / pixel, what dilution of beads do you need to make to achieve the desired density in the viewing area? The density of polystyrene is ≈ 1.05 g/ml. Assume that the depth of field ≈ 5 μm and that you are doing your first measurements with beads with diameter 1 μm . Assume, also, that the concentration is as specified in Sec. 2.2.

Pixel pitch:

3.45 μm per pixel

WRONG - 100nm per pixel

Conc. 0.5% weight by volume

20 beads per 1440x1080 pixels x 0.1 μm per pixel
= 20 beads per 144 x 108 μm x 5 μm

G of beads

0.5 μm rad per bead = 0.5236 μm^3

/1000³ = 0.5236 * 10⁻⁹ mm³

/1000 = 0.5236 * 10⁻¹² cm³ (ml)

X 1.05 g/ml = 5.4978 * 10⁻¹³ g/bead

Volume of image = 77,760 μm^3 / 1000³ = 77.76 * 10⁻⁶ mm³

/1000 = 7.776 * 10⁻⁸ ml

Beads per ml = 20 / 7.776 * 10⁻⁸ ml = 2.572 * 10⁸ beads/ml

g/bead * bead/ml = g/ml = 0.0001414

0.005g/1ml initial

Initial/final conc. = 35.36X dilution
