

Lab 2 Session 1

February 3, 2026 1:06 PM

Notes

- Fix assignment 1
- Everything on Goodnote?
- Snapshot 1 due Thursday
- Assignment 2 due feb 10, 3 due feb 24 after break

2026-02-03 1:31 PM

Lab 2 Microscopy and Motility

Nathan Unruh

Ahilan

Table

- Only collaborations in lab not outside (lab notes can be done together but reflections data ect is individual)

Physical image length per pixel calibration:

8 lines, 1168.706

73.044 pixels per line

From slide: 0.01mm per line

$\text{Mm/px} = 6.845177 \times 10^{-5} \text{mm} \times 10^3 = 0.06845177 \text{ um / px} = 68.45177 \text{nm/px}$

Camera model BFS-U3-16S2C-CS

Calculated dilution factor (using prelab method)

20 beads per 1440x1080 pixels x 0.0685um per pixel

= 20 beads per 98.6 x 73.98um x 5um

Volume of image: 36,482 um^3 vs 77,760 um^3 = 0.46916x
(compared to prelab)

Density Prelab: 2.572×10^8 beads/ml

Density our lab: 1.207×10^8 beads/ml

5.4978×10^{-13} g/bead * 1.207×10^8 beads/ml = 6.636×10^{-5} g/ml

Initial: 0.005g/ml

Dilution = initial/final = 75.34X